

## Pain Management Comparison Table



Model	E4	E3 Intense	E2 ELITE
<b>Product code</b>	HV-F128-E	HV-F021-EW	HV-F127-E
<b>Triple Action TENS<sup>5</sup> Technology</b> Blocks the pain message, triggers the release of endorphins, improves the blood circulation	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input checked="" type="radio"/>
<b>Drug-Free</b>	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input checked="" type="radio"/>
<b>Portable - convenient to take with you</b>	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>
<b>Treatment zones</b>	6: back, shoulder, joints, leg, arm, sole (or feet)	6: back, shoulder, joints, leg, arm, sole (or feet)	4: back, shoulder, leg, sole (or feet)
<b>Massage modes</b>	4: tap, knead, rub, push	3: tap, knead, rub	4: tap, knead, rub, push
<b>Intensity levels</b>	10	15	10
<b>Quick Relief button</b> Uses high frequency vibrations to relieve acute pain promptly	<input checked="" type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>
<b>Hyper mode</b> Increases from gradual to a higher strength, helping relief pain effectively	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
<b>Soft mode</b> Gentle program which helps prevent against aggravating existing pain	<input checked="" type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>
<b>Selection of left/right pad treatment</b>	<input checked="" type="radio"/>	<input type="radio"/>	<input checked="" type="radio"/>
<b>Display</b>	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>
<b>Pads</b>	Long Life pads, reusable up to 150 times	Long Life pads, reusable up to 150 times	Long Life pads, reusable up to 150 times

## Pain Relievers

### Accessories

Model	Product code	Description
Long Life Pads	HV-LLAPD-E	Whashable pads, reusable up to 150 times Can be used with OMRON Pain Relievers E4, E3 Intense and E2 Elite