

Weight Management Comparison Table



Model	BF511	BF508	BF214	BF212	BF306	HN288	HN286	HN289
Product Code	HBF-511T-E HBF-511B-E	HBF-508-E	HBF-214-EBW	HBF-212-EW	HBF-306-E	HN-288-E	HN-286-E	HN-289-ESL HN-289-EB HN-289-EPK HN-289-EBK
Body Fat (in %) Amount of body fat mass as percentage of total body weight	⊙	⊙	⊙	⊙	⊙	○	○	○
Body Mass Index Ratio between the weight and height of a person	⊙	⊙	⊙	⊙	⊙	○	○	○
Skeletal Muscle (in %) Muscles attached to bones that move the body	⊙	○	⊙	○	○	○	○	○
Visceral Fat Fat surrounding the organs	⊙	⊙	○	○	○	○	○	○
Resting Metabolism Rate Minimum level of caloric intake needed to sustain the body's everyday functions	⊙	○	○	○	○	○	○	○
Full Body Measurement Both hand and foot measurement with 8 sensors	⊙	⊙	○	○	○	○	○	○
Clinically Validated⁷	⊙	⊙	⊙	⊙	⊙	○	○	○
Number of user profiles that can be set up (in addition to guests)	4	4	4	4	9	○	○	○
Memory function Shows previous measurement	⊙	○	⊙	○	○	○	○	○
Weight difference Effortless weigh smaller things like a baby or luggage	○	○	○	○	○	⊙	○	○
Age range	6-80 yrs	18-80 years	10-80 years	10-80 years	10-70 years	n/a	n/a	n/a
Increments	100 gr	100 gr	100 gr	100 gr	100 gr	100 gr	100 gr	100 gr
Maximum weight capacity	150 kg	150 kg	150 kg	150 kg	199.8 kg	180 kg	180 kg	150 kg
Material	stainless steel/ plastic	stainless steel/ plastic	stainless steel/ glass	stainless steel/ glass	stainless steel/ plastic	glass/ plastic	glass/ plastic	glass/ plastic
On/Off technology	switch	switch	sensor switch	sensor switch	switch	auto on/off	auto on/off	auto on/off